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Equipment on Counter-Balanced Arms

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One Time Checks (12 months after installation)

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WORK DETAILS

Engineer Name: TAKATOMO SAKAI

Event Type: Predictive / Preventative Maintenance

Device Status: Pass

Service Activity: Visit 1

EQUIPMENT DETAILS

Product Name: Azurion 7 B12

Serial Number: 19

Product Number: 722225

Installed Product: 88127914

UDI: (01)00884838099265(21)19

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

On-site preparations

Preparing on-site for Planned Maintenance

Schedule a patient	PASS
Check emergency stop	PASS
Check safety of system	PASS
Check system	PASS

Image Chain + X-ray Safety

Image chain Frontal channel

Adjusting the generator	PASS
Check Dose Test and Fastdose Calibration	PASS
Check Energy and Signal Flow 1	PASS
Check Energy and Signal Flow 2	PASS
Check EDL Expected Tube Power Verification	PASS

X-ray safety: Air Kerma rate Verification (frontal)

Dose Rate Fluoroscopy Mode I	μGy/s	374	
Local limit value: Fluoroscopy Mode I			833 μGy/s
Dose Rate Fluoroscopy Mode I within local limit			PASS
Dose Rate Fluoroscopy Mode II	μGy/s	748	
Local limit value: Fluoroscopy Mode II			833 μGy/s
Dose Rate Fluoroscopy Mode II within local limit			PASS
Dose Rate Fluoroscopy Mode III	μGy/s	1882	
Local limit value: Fluoroscopy Mode III			2083 μGy/s

	Min	Max	UoM	Measurement	Result	Comments
Image Chain + X-ray Safety						
X-ray safety: Air Kerma rate Verification (frontal)						
Dose Rate Fluoroscopy Mode III within local limit					PASS	
Image chain Lateral channel						
Adjusting the generator					PASS	
Check Dose Test and Fastdose Calibration					PASS	
Check Energy and Signal Flow 1					PASS	
Check Energy and Signal Flow 2					PASS	
Check EDL Expected Tube Power Verification					PASS	
X-ray safety: Air Kerma rate Verification (lateral)						
Dose Rate Fluoroscopy Mode I			μGy/s	370		
Local limit value: Fluoroscopy Mode I					833 μGy/s	
Dose Rate Fluoroscopy Mode I within local limit					PASS	
Dose Rate Fluoroscopy Mode II			μGy/s	738		
Local limit value: Fluoroscopy Mode II					833 μGy/s	
Dose Rate Fluoroscopy Mode II within local limit					PASS	
Dose Rate Fluoroscopy Mode III			μGy/s	1860		
Local limit value: Fluoroscopy Mode III				2083 μGy/s		
Dose Rate Fluoroscopy Mode III within local limit					PASS	
Image chain Monitor (Philips)						
Check Monitor Performance					PASS	
Mechanical Maintenance + Preventative Replacements						
Operator controls						
Check the function of the operator controls					PASS	
Calibrations						
Interventional Tools / XtraVision: XperCT gain calibration						
Do the XperCT gain calibration					Not Required	
Completing the visit						
Completing the visit						
Resetting system logging					PASS	
Handover					PASS	

PHILIPS

Philips Representative Signature:

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Engineer Name: TAKATOMO SAKAI

Date: Aug 30, 2021

Page: 3 / 3

Work Order: WO-04054584