

保守点検報告書

2021年 第1回 保証期間内点検

顧客名 長崎大学病院
装置名 Azurion7B12
設置場所 血管造影室31室
点検実施日 2021年8月23日
Case番号 0118005845

株式会社フィリップス・ジャパン

Ship To:
 長崎大学病院
 坂本1-7-1
 長崎市
 Japan

Test and Verification Report

Work Order: WO-04054582

WORK DETAILS

Engineer Name: TAKATOMO SAKAI
 Device Status: Pass

Event Type: Predictive / Preventative Maintenance
 Service Activity: Visit 1

EQUIPMENT DETAILS

Product Name: Azurion 7 B12
 Product Number: 722225

Serial Number: 18
 Installed Product: 88127912
 UDI: (01)00884838099265(21)18

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

On-site preparations

Preparing on-site for Planned Maintenance

Schedule a patient					PASS	
Check emergency stop					PASS	
Check safety of system					PASS	
Check system					PASS	

Image Chain + X-ray Safety

Image chain Frontal channel

Adjusting the generator					PASS	
Calibrate the flat detector					PASS	
Check Dose Test and Fastdose Calibration					PASS	
Check Energy and Signal Flow 1					PASS	
Check Energy and Signal Flow 2					PASS	
Check EDL Expected Tube Power Verification					PASS	

X-ray safety: Air Kerma rate Verification (frontal)

Dose Rate Fluoroscopy Mode I			μGy/s	370		
Local limit value: Fluoroscopy Mode I					833 μGy/s	
Dose Rate Fluoroscopy Mode I within local limit					PASS	
Dose Rate Fluoroscopy Mode II			μGy/s	737		
Local limit value: Fluoroscopy Mode II					833 μGy/s	
Dose Rate Fluoroscopy Mode II within local limit					PASS	
Dose Rate Fluoroscopy Mode III			μGy/s	1855		

	Min	Max	UoM	Measurement	Result	Comments
Image Chain + X-ray Safety						
X-ray safety: Air Kerma rate Verification (frontal)						
Local limit value: Fluoroscopy Mode III				2083 $\mu\text{Gy/s}$		
Dose Rate Fluoroscopy Mode III within local limit					PASS	
Image chain Lateral channel						
Adjusting the generator					PASS	
Calibrate the flat detector					PASS	
Check Dose Test and Fastdose Calibration					PASS	
Check Energy and Signal Flow 1					PASS	
Check Energy and Signal Flow 2					PASS	
Check EDL Expected Tube Power Verification					PASS	
X-ray safety: Air Kerma rate Verification (lateral)						
Dose Rate Fluoroscopy Mode I			$\mu\text{Gy/s}$	369		
Local limit value: Fluoroscopy Mode I					833 $\mu\text{Gy/s}$	
Dose Rate Fluoroscopy Mode I within local limit					PASS	
Dose Rate Fluoroscopy Mode II			$\mu\text{Gy/s}$	738		
Local limit value: Fluoroscopy Mode II					833 $\mu\text{Gy/s}$	
Dose Rate Fluoroscopy Mode II within local limit					PASS	
Dose Rate Fluoroscopy Mode III			$\mu\text{Gy/s}$	1860		
Local limit value: Fluoroscopy Mode III				2083 $\mu\text{Gy/s}$		
Dose Rate Fluoroscopy Mode III within local limit					PASS	
Image chain Monitor (Philips)						
Check Monitor Performance					PASS	
Mechanical Maintenance + Preventative Replacements						
Frontal stand						
Do the mechanical maintenance of the frontal stand					PASS	
Operator controls						
Check the function of the operator controls					PASS	
Calibrations						
Frontal stand						
Adjust frontal stand					PASS	
Interventional Tools / XtraVision: XperCT gain calibration						
Do the XperCT gain calibration					Not Required	

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

Calibrations

Interventional Tools / XtraVision: calibrations and verifications

Do the Interventional Tools / XtraVision calibrations and verifications

PASS

Completing the visit

Completing the visit

Resetting system logging

PASS

Handover

PASS

Philips Representative Signature:

酒井 崇友

Engineer Name: TAKATOMO SAKAI

Date: Aug 23, 2021

